




 **SYRINGES**  
LARGE CAPACITIES





**PENTA “Large Capacities”** syringes are dedicated to those special needs that require a secure grip and a high level of precision even in the presence of large amounts of liquid to be injected. The 50/60 ml syringe is also available in opaque amber polypropylene to ensure protection for photosensitive medication. The 100 ml syringes are fitted with a Luer adapter, which is detachable and located in the plunger. They have **“Latex Free”** piston, made of synthetic rubber which does not contain natural latex, in order to prevent any allergies. Two types: with and without needle.

### Three different types of cone

- The Eccentric Luer cone facilitates the operator during sampling.
- The Luer Lock cone makes the connection with the needle and/or additional connector ensuring a tight seal.
- The catheter cone enables the attachment of probes.



Adapter 100 ml

### Characteristics of the components

- Polypropylene **cylinder** with retaining ring.
- Polypropylene **plunger**.
- **“Latex Free” piston** with two seal rings, made of non-toxic, latex-free synthetic rubber.
- AISI 304 steel **cannula**; diameter in accordance with ISO 9626 and length in accordance with ISO 7864.
- Coloured polypropylene **hub** in accordance with ISO 6009.
- Transparent polypropylene **100ml Luer Adapter**.
- Transparent polypropylene **needle cover**.
- Medical grade silicone **lubricant**.
- Indelible **graduation** in accordance with ISO 7886-1.
- Medical grade **blister** paper and plastic film.

### Sterilization method

Ethylene Oxide (EO 10% + CO2 90%), cycle validated in accordance with ISO 11135-1

#### Ferrara site

Via Modena 119,  
44122 Ferrara  
Tel. +39.0532.778811  
Fax +39.0532.778812

#### Campli site

S.S. 202 Loc.tà Nocella,  
64012 Campli (TE)  
Tel. +39.0861.560201  
Fax +39.0861.560200

info@pentaferite.com  
www.pentaferite.com

CAPACITY	GRADUATION	NOZZLE	NEEDLE	PACKAGING
30 ml	1 ml	B - C	L - S	80 pcs
50 ml *	1 ml	C	L - S	60 pcs
50/60 ml	1 ml	B - C	L - S	60 pcs
50/60 ml	1 ml	D	S	55 pcs
100 ml	2 ml	D **	-	50 pcs

A: Central Luer  
B: Eccentric Luer  
C: Luer Lock  
D: Catheter

L: With needle on side  
S: Without needle  
\* Amber  
\*\* With Luer adapter