

ACCESSORIES
FOR ENTERAL FEEDING



ACCESSORIES

FOR ENTERAL FEEDING



Pentaferte has a wide range of accessories for enteral feeding which, along with the Penta Segap syringes, complete a range of products which enable solutes and medication to be prepared and administered in complete safety.

Accessory features

Medicine straw: device to draw the medication from the medicament bottle

- Purple male **connector** or magenta WD.
- Medical grade DEHP Free soft PVC **tube**:
 - inner diameter 1 mm
 - outer diameter 2 mm
 - length 5 cm
- Medical grade **blister** paper and plastic film.

Milk straw: device to draw milk or solutes from the bottle

- **Purple** male connector.
- Medical grade transparent soft PVC **tube**:
 - inner diameter 2.9 mm
 - outer diameter 4 mm
 - length 20 cm
- Medical grade **blister** paper and plastic film.

Bolus filter: device for closing 20/60 ml syringes to prevent infection

- Purple **cap** with disc filter.
- Medical grade **blister** paper and plastic film.

3-way stopcock: made of polycarbonate/polyethylene plastic material

- **Connector:** 2 male Luer Lock, 1 female Luer Lock.
- Medical grade **blister** paper and plastic film.

Purple **male Luer Lock cap**, can be adapted to any Penta Segap syringe with female Luer cone

- Medical grade **blister** paper and plastic film.

Nasogastric tubes with low residue and male connector

- Purple male **connector**.
- **Tube** available in DEHP Free PVC, soft polyurethane or medical grade transparent silicone:
 - graduation every 1 cm
 - two side holes opposite each other; positioned and shaped in an asymmetrical way to prevent the risk of occlusion
 - radio-opaque band for X-rays
 - available in sizes CH 4-5-6-8-10-12
- Coloured cap-holder **tab**.
- Connector closing **cap**.
- Medical grade **blister** paper and plastic film.

Extension tube: made of non-toxic medical grade DEHP Free PVC.

- Purple **male connector** or magenta WD connector.
- Purple **female connector** or magenta WD connector.
- Medical grade transparent soft PVC **tube** with low residual volume:
 - inner diameter 1.40 mm
 - outer diameter 2.67 mm
 - length 150 cm / 90 cm
 - residual volume 2.30 ml (x 150) / 1.28 ml (x 90)
- Medical grade **blister** paper and plastic film.

Sterilization method

Ethylene Oxide (EO 10% + CO₂ 90%), cycle validated in accordance with ISO 11135-1

Ferrara site

Via Modena 119,
44122 Ferrara
Tel. +39.0532.778811
Fax +39.0532.778812

Campli site

S.S. 202 Loc.tà Nocella,
64012 Campli (TE)
Tel. +39.0861.560201
Fax +39.0861.560200

info@pentaferte.com
www.pentaferte.com