

PentaFerte

MEDICAL DEVICES SOLUTIONS



1ML POLYCARBONATE LUER LOCK SYRINGE

Features

- ✓ **STERILE SINGLE USE MEDICAL DEVICE**
- ✓ **DEHP FREE / LATEX FREE / PVC FREE**
- ✓ **GRADUATION SCALE INTERVAL: 0,01 ML**
- ✓ **LUER LOCK TIP**
- ✓ **DOUBLE SEALING RING GASKET**
- ✓ **ETHYLENE OXIDE STERILIZATION**



Polycarbonate resin provides **reliable connections** avoiding the risk of leakage or needle detachment

Barrel clearness facilitates
fluid monitoring



Additional info

- ✓ **SHELF LIFE: 5 YEARS**
- ✓ **CASE: 1.800 UNITS / 18 BOXES**
- ✓ **BOX: 100 UNITS**
- ✓ **REF: 002022120P**