

PentaFerte

MEDICAL DEVICES SOLUTIONS



LARGE CAPACITY SYRINGES

About Large Capacity Syringes

Penta large capacity range meets **safety and accuracy** needs even in case of **high volumes to be injected**.

50/60 ml syringe is also available with amber light resistant version to guarantee light sensitive drugs protection.

100 ml capacity is supplied with detachable Luer adapter housed in the plunger.

Three different nozzles:

- ✓ **Eccentric Luer nozzle** facilitates the operator with samples drawing
- ✓ **Luer Lock nozzle** provides tight reliable connections
- ✓ **Catheter nozzle** allows tubes connection

Components specifications

- ✓ **Barrel:** transparent polypropylene with plunger retention ring
- ✓ **Plunger:** opaque polypropylene
- ✓ **Plunger stopper:** synthetic non-toxic rubber with two sealing rings
- ✓ **Luer adapter** (100 ml capacity only): transparent polypropylene
- ✓ **Lubricant:** medical grade silicone
- ✓ **Graduation scale:** according to EN ISO 7886-1
- ✓ **Unit packaging (sterile barrier system):** medical grade paper - plastic film

Sterilization method

Ethylene Oxide 10% + CO2 90%; cycle validated according to EN ISO 11135.

| Code | Description | Nozzle | Units/Carton | Boxes/Carton |
|-------------------------|------------------------------------|--------|--------------|--------------|
| 002022810 | SYR.30ML ECC. | B | 320 | 4 X 80 PCS |
| 002022820F | SYR.30ML L/L | C | 320 | 4 X 80 PCS |
| 002022920AF | SYR.50ML L/L AMBER-LIGHT PROTECTED | C | 240 | 4 X 60 PCS |
| 002022960 | SYR.60ML ECC. | B | 240 | 4 X 60 PCS |
| 002022970 002022970F | SYR.60ML L/L | C | 240 | 4 X 60 PCS |
| 002022980 | SYR.60ML CATHETER | D | 220 | 4 X 55 PCS |
| 002022C30 | SYR.100ML CATH. W/LUER ADAPTOR* | D | 50 | 50 X 1 PCS |

Additional info

- B:** Eccentric Luer
- C:** Luer Lock
- D:** Catheter

