

PentaFerte

MEDICAL DEVICES SOLUTIONS



BLOOD LANCETS

About Blood Lancets

Designed to facilitate blood draw. The ultra sharp tip allows an easy penetration while the plastic handle allows easy handling of the device.

General specifications

- ✓ Sterile, single use, non-toxic

Sterilization method

Gamma rays 25 KGy; cycle validated according to EN ISO 11137

Components specifications

- ✓ **Lancet handle:** polyethylene length 23mm (+/-0.1mm), width 5mm (+/-0.1mm)
- ✓ **Tip:** 30G, stainless steel AISI 304, primary angle of 11°+/-2° according to EN ISO 7864, ext. 0.30mm (according to EN ISO 9626)
- ✓ **Tip protector:** polyethylene 9.4mm (+/-0.1mm)
- ✓ **Packaging:** white carton box

Code	Description	Units/Carton	Boxes/Carton
1290210301	Blood Lancets	20.000	200 x 100 PCS

