



## General specifications

The Penta Enteral accessories range with WD connector completes the range of Penta Enteral Syringes and is designed for the preparation and administration of enteral fluids and drugs to the gastrointestinal system of newborns, pediatric patients and adults. The WD connector, compliant with EN ISO 80369-1 and EN ISO 1615 standards, reduces the risk of accidental connection to a parenteral system with Luer or Luer Lock adapter. The purple colour allows the differentiation from lines of accessories and syringes with different intended use.

## Sterilization method

Ethylene Oxide 10% + CO<sub>2</sub> 90%; cycle validated according to EN ISO 11135.

## Packaging:

Unit packaging (sterile barrier system): medical grade paper - plastic film.

Code	Description	Units / Carton	Boxes / Carton
032022AC02WD	MEDICINE STRAW 5CM WD PENTAENTERAL	200	N/A
Straw for medicine drawing up from container: purple colour-coded WD connector; Soft medical grade PVC; Residual volume 0,10 ml; Length 5 cm.			
032022AC01WD	MILK STRAW 20CM WD PENTAENTERAL	200	N/A
Straw for milk or nutritional formula drawing up from container: purple colour-coded WD connector; Soft medical grade PVC; Residual volume 1,30 ml; Length 20 cm.			
032022AC08WD	BOLUS FILTER CAP FOR SYRINGES WD PENTAENTERAL	150	N/A
Purple colour-coded bolus filter cap 15 micron to reduce the risk of bacterial infection for 20/60 ml syringes. Male ended 3-way stopcock allows line flushing and drug delivery ; 2 WD Male connectors, 1 WD Female connector.			
032022AC03WD	MALE CAP L/L 1 BLISTER OF 1 WD PENTAENTERAL	200	N/A
032022AC04WD	MALE CAP L/L 1 BLISTER OF 8 WD PENTAENTERAL	200	N/A
Purple colour-coded WD male cap compatible with any Penta syringe with female WD connector.			
032752SC04WD	NASO-GASTRIC TUBE PVC CH5/50 WD PENTAENTERAL	400	8 X 50 PZ
032752SC06WD	NASO-GASTRIC TUBE PVC CH6/50 WD PENTAENTERAL	400	8 X 50 PZ
032752SC06WD1	NASO-GASTRIC TUBE PVC CH6/80 WD PENTAENTERAL	400	8 X 50 PZ
032752SC06WD2	NASO-GASTRIC TUBE PVC CH6/120 WD PENTAENTERAL	400	8 X 50 PZ
032752SC08WD	NASO-GASTRIC TUBE PVC CH8/50 WD PENTAENTERAL	400	8 X 50 PZ
032752SC08WD1	NASO-GASTRIC TUBE PVC CH8/80 WD PENTAENTERAL	400	8 X 50 PZ
032752SC08WD2	NASO-GASTRIC TUBE PVC CH8/120 WD PENTAENTERAL	400	8 X 50 PZ
032752SC10WD2	NASO-GASTRIC TUBE PVC CH10/120 WD PENTAENTERAL	400	8 X 50 PZ
032752SH04WD	NASO-GASTRIC TUBE PUR CH5/50 WD PENTAENTERAL	400	8 X 50 PZ
032752SH06WD	NASO-GASTRIC TUBE PUR CH6/50 WD PENTAENTERAL	400	8 X 50 PZ
032752SH06WD3	NASO-GASTRIC TUBE PUR CH6/120 WD PENTAENTERAL	400	8 X 50 PZ
032752SH08WD	NASO-GASTRIC TUBE PUR CH8/50 WD PENTAENTERAL	400	8 X 50 PZ
032752SH08WD3	NASO-GASTRIC TUBE PUR CH8/120 WD PENTAENTERAL	400	8 X 50 PZ
032752SH10WD2	NASO-GASTRIC TUBE PUR CH10/120 WD PENTAENTERAL	320	8 X 40 PZ
Low residual Naso-gastric tubes: purple colour-coded Male connector; Available in PVC and medical grade PUR; Graduation every 1 cm; Two opposite lateral holes, positioned and shaped asymmetrically to avoid the risk of occlusion; X-ray radio-opaque stripe; Available sizes: CH 4-5-6-8-10-12; Colored cap-holder; Connector closing cap.			
032022AC06WD2	EXTENSION SET 90CM WD PENTAENTERAL	200	N/A
032022AC07WD	EXTENSION SET 150CM WD PENTAENTERAL	200	N/A
Extension sets: Purple colour-coded Male and Female WD connectors; Soft medical grade PVC tube; Low residual volume. ID: 1,40 mm: Residual volume 2,30 ml (150 cm length) / 1,38 ml (90 cm length).			
032022AC10WDC	ADAPTER ADRS-CATH WD PENTAENTERAL	200	N/A
032022AC11WD	ADAPTER FOR NASO-GASTRIC TUBE ADRS WD PENTAENTERAL	200	N/A
032022AC13WD	SYRINGE BOTTLE ADAPTER MEDIUM WD PENTAENTERAL	150	N/A