

PentaFerte

MEDICAL DEVICES SOLUTIONS



PENTAFINE

About Three-Piece R0 Syringes

PentaFine R0 (minimized Dead Space) range represents an immediate and non-traumatic solution to address diabetic patients' needs ensuring total and painless injection.

General specifications

- ✓ **Sterile**, single use, non-toxic
- ✓ **Dead space** in compliance with EN ISO 8537 specifications
- ✓ **Fixed needle** ensures safe painless injection and accurate administration
- ✓ **Pad-printed graduation** scale allows easier reading
- ✓ **Available capacities:**
0,3ml 100-UI, 0,5ml 100-UI.

Sterilization method

Ethylene Oxide.

Components specifications

- ✓ **Barrel:** transparent polypropylene
- ✓ **Plunger:** opaque polypropylene
- ✓ **Plunger stopper:** synthetic non-toxic rubber with two sealing rings
- ✓ **Needle cannula:** stainless steel SUS 304; diameter and length according to EN ISO 7864.
- ✓ **Needle cap:** polypropylene plus orange colourant
- ✓ **Lubricant:** medical grade silicone
- ✓ **Graduation scale:** according to EN ISO 8537
- ✓ **Unit packaging (sterile barrier system):** medical grade paper - plastic film

Code	Description	Needle	Units/Carton	Boxes/Carton
00203202P03	SYR. 0,3ML 100-UI R.O.	30G x 1/3" 0,30x8mm	2.400	24 X 100 PCS
00203207N05	SYR. 0,5ML 100-UI R.O.	29G x 1/2" 0,33x12mm	2.400	24 X 100 PCS
00203207P03	SYR. 0,5ML 100-UI R.O.	30G x 1/3" 0,30x 8mm	2.400	24 X 100 PCS
00203207P05	SYR. 0,5ML 100-UI R.O.	30G x 1/2" 0,30x12mm	2.400	24 X 100 PCS