

PentaFerte

MEDICAL DEVICES SOLUTIONS



THREE-PIECE SYRINGES WO/NEEDLE

About Three-Piece Syringes wo/needle

Penta syringes offer a complete range of products manufactured in compliance with high-quality and safety standards.

General specifications

- ✓ **Sterile**, single use, non-toxic
- ✓ Barrel wall material enables **easy dosage monitoring** and greater control over mixing and liquid flow
- ✓ **Black graduation scale** for easier reading
- ✓ Medical grade silicone oil treatment guarantees **smooth and safe piston sliding**
- ✓ Push-button/Finger grips distance according to EN ISO 7886-1 guarantees **safe and easy handling**
- ✓ **Luer and Luer Lock nozzles** designed according to EN ISO 80369-7

Components specifications

- ✓ **Barrel:** transparent polypropylene with plunger retention ring
- ✓ **Plunger:** opaque polypropylene
- ✓ **Plunger stopper:** synthetic non-toxic rubber with two sealing rings
- ✓ **Lubricant:** medical grade silicone
- ✓ **Graduation scale:** according to EN ISO 7886-1
- ✓ **Unit packaging (sterile barrier system):** medical grade paper - plastic film

Sterilization method

Ethylene Oxide 10% + CO2 90%; cycle validated according to EN ISO 11135.

Code	Description	Nozzle	Units/Carton	Boxes/Carton
002022120	SYR.1ML INS.100-UI	A	2.400	24 X 100 PCS
002022140	SYR.1ML TUB.	A	2.400	24 X 100 PCS
002022142	SYR.1ML L/L	C	2.000	20 X 100 PCS
002022300	SYR.2,5ML CEN.	A	2.000	20 X 100 PCS
002022320	SYR.2,5ML L/L	C	2.000	20 X 100 PCS
002022500	SYR.5ML CEN.	A	1.600	16 X 100 PCS
002022520F	SYR.5ML L/L	C	1.600	16 X 100 PCS
002022600	SYR.10ML CEN.	A	1.000	10 X 100 PCS
002022610	SYR.10ML ECC.	B	1.000	10 X 100 PCS
002022620F	SYR.10ML L/L	C	1.000	10 X 100 PCS
002022710	SYR.20ML ECC.	B	800	10 x 80 PCS
002022720F	SYR.20ML L/L	C	800	10 x 80 PCS
002022720AF	SYR.20ML L/L AMBER-LIGHT PROTECTED*	C	800	10 x 80 PCS

Additional info

- A: Central Luer
- B: Eccentric Luer
- C: Luer Lock



PentaFerte Italia Srl Unipersonale - VAT: IT01866900671
Registered Office & Production Site: Viale Piane Nocella, 23 - 64012 Campli (TE) - Italy
Commercial & Administrative Offices: Via Modena, 119/A - 44122 Ferrara - Italy