

PentaFerte

MEDICAL DEVICES SOLUTIONS



THREE-PIECE SYRINGES W/NEEDLE

About Three-Piece Syringes w/needle

Penta syringes offer a complete range of products manufactured in compliance with high-quality and safety standards.

General specifications

- ✓ Sterile, single use, non-toxic
- ✓ Barrel wall material enables **easy dosage monitoring** and greater control over mixing and liquid flow
- ✓ **Black graduation scale** for easier reading
- ✓ Medical grade silicone oil treatment guarantees **smooth and safe piston sliding**
- ✓ Push-button/Finger grips distance according to EN ISO 7886-1 guarantees **safe and easy handling**
- ✓ **Luer nozzle** designed according to EN ISO 80369-7

Sterilization method

Ethylene Oxide 10% + CO₂ 90%; cycle validated according to EN ISO 11135.

Components specifications

- ✓ **Barrel:** transparent polypropylene with plunger retention ring
- ✓ **Plunger:** opaque polypropylene
- ✓ **Plunger stopper:** synthetic non-toxic rubber with two sealing rings
- ✓ **Needle cannula:** stainless steel AISI 304; diameter according to EN ISO 9626 and length according to EN ISO 7864
- ✓ **Needle hub:** coloured polypropylene according to EN ISO 6009
- ✓ **Needle cap:** transparent polypropylene
- ✓ **Lubricant:** medical grade silicone
- ✓ **Graduation scale:** according to EN ISO 7886-1
- ✓ **Unit packaging (sterile barrier system):** medical grade paper - plastic film

Code	Description	Nozzle	Needle	Units/Carton	Boxes/Carton
00202212H58	SYR.1ml INS.100U 25G x 5/8"	A	0,50 X 16mm	2.000	20 X 100 PCS
00202212I05	SYR.1ml INS.100U 26G x 1/2"	A	0,45 X 12mm	2.000	20 X 100 PCS
00202212L05	SYR.1ml INS.100U 27G x 1/2"	A	0,40 X 12mm	2.000	20 X 100 PCS
00202212M05	SYR.1ml INS.100U 28G x 1/2"	A	0,36 X 12mm	2.000	20 X 100 PCS
00202213H58U	SYR.1ml INS/TUB 100U 25G x 5/8"	A	0,50 X 16mm	2.000	20 X 100 PCS
00202213I05U	SYR 1ml INS/TUB 100U 26G x 1/2"	A	0,45 X 12mm	2.000	20 X 100 PCS
00202214H58	SYR.1ml TUB. 25G x 5/8"	A	0,50 X 16mm	2.000	20 X 100 PCS
00202214I05	SYR.1ml TUB. 26G x 1/2"	A	0,45 X 12mm	2.000	20 X 100 PCS
00202230D12	SYR.2,5ml 21G x 1 1/2"	A	0,80 X 40mm	2.000	20 X 100 PCS
00202230E14	SYR.2,5ml 22G x 1 1/4"	A	0,70 X 30mm	2.000	20 X 100 PCS
00202230F10	SYR.2,5ml 23G x 1"	A	0,60 X 25mm	2.000	20 X 100 PCS
00202230F14	SYR.2,5ml 23G x 1 1/4"	A	0,60 X 30mm	2.000	20 X 100 PCS
00202250D12	SYR.5ml CEN. 21G x 1 1/2"	A	0,80 X 40mm	1.200	12 X 100 PCS
00202250E12	SYR.5ml CEN. 22G x 1 1/2"	A	0,70 X 40mm	1.200	12 X 100 PCS
00202250E14	SYR.5ml CEN. 22G x 1 1/4"	A	0,70 X 30mm	1.200	12 X 100 PCS
00202250F14	SYR.5ml CEN. 23G x 1 1/4"	A	0,60 X 30mm	1.200	12 X 100 PCS
00202260E12	SYR.10ml CEN. 22G x 1 1/2"	A	0,70 X 40mm	960	12 X 80 PCS
00202260E14	SYR.10ml CEN. 22G x 1 1/4"	A	0,70 X 30mm	960	12 X 80 PCS
00202260D12	SYR.10ml CEN. 21G x 1 1/2"	A	0,80 X 40mm	960	12 X 80 PCS
00202261D12	SYR.10ml ECC. 21G x 1 1/2"	B	0,80 X 40mm	960	12 X 80 PCS
00202271C12L	SYR.20ml ECC. 20G x 1 1/2" (with side needle)	B	0,90 X 40mm	800	10 X 80 PCS
00202271D12L	SYR.20ml ECC. 21G x 1 1/2" (with side needle)	B	0,80 X 40mm	800	10 X 80 PCS
00202271E14L	SYR.20ml ECC. 22G x 1 1/4" (with side needle)	B	0,70 X 30mm	800	10 X 80 PCS

Additional info

Non-Standard
Combination Syringes/
Needle: available upon
request

A: Central Luer
B: Eccentric Luer